



Hahnemühle

“HaMuNa® Care” – FFP2 Mask Made in Germany

Flu waves and the corona pandemic have shown us how important personal protection against infections with viruses and bacteria is. A filtering half mask is the best every day protection in all public places.

The newly developed half mask from Hahnemühle with excellent filter performance represents personal protective equipment (PPE) and protects the carrier from particles, droplets, and aerosols. The air becomes filtered by both in and exhalation and thus also serves as external protection. Is comfortable to wear, non irritating and seals very well. Fogging of glasses is thereby significantly reduced. The material is odour-free with a pleasant surface.



Components are:

- Three high particle filtering fleece layers, thermally bonded polypropylene fibres
- Individually adjustable galvanized nosepiece with HDPE coating
- Elastic bands made of polyamide and elastane

Features:

- Approval as PPE according to DIN EN149: 2001 + A1: 2009 –FFP2 NR
- Almost 98% filter efficiency against solid and liquid aerosols for optimal protection
- Additional quality label “Covid Certified Filter (CCF) 0158 150 Nano”
- CE marking

Characteristics:

- Light, high-quality fabric for comfort wearing
- Individually adjustable nose clip for greatly reduced fogging of the glasses
- Elastic, comfortable rubber bands
- Odour-free materials

Important: The mask does not protect against gases and vapours.

Durability:

2 years from the date of manufacture when stored in the original packaging, away from contaminated areas, dust, sunlight, temperature extremes, excessive moisture, and harmful chemicals.

Recommended storage temperatures between - 20 ° C and + 40 ° C and relative humidity of less than 60%.

All HaMuNa® Care products are characterized by the fact that they are exclusively manufactured in Germany. Therefore, we can guarantee a high-quality standard especially in this sensitive area.

As a global company, Hahnemühle focuses on products in the areas of life science, including healthcare, as well as chemical analysis. Together with customers, in a relentless pursuit of innovation, it continues to develop new products that are traded at the highest level on the world market. With production processes that are constantly being further developed and strategic partnerships with selective expertise, we serve our customers with new products, the reproducible quality of which gives our customers a competitive edge.